

Carestream

DRX-Revolution Nano Mobile X-ray System

Small. Agile. Affordable.

Revolutionary Performance on a Smaller Scale.



SMART IMAGING SOLUTIONS

Striking the Ideal Balance of Price and Performance.

Looking for a mobile X-ray system to replace your existing unit or add to your current fleet? CR systems, while affordable, don't provide all the benefits of a fully digital system. At the same time, adding additional, fully integrated digital systems may be out of your price range.

Meet the best of both worlds. We've leveraged the advanced features and digital integration of our flagship DRX-Revolution to create the small and compact DRX-Revolution Nano Mobile X-ray System. It delivers high performance and DR imaging, scaled down to make it an economical choice. It's just the Revolution you've been waiting for.

Advantages at a Glance

- Small and lightweight for easy transport – even in tight, cluttered environments.
- Balanced, articulating arm for extended reach and easy positioning.
- Carbon Nano Tube to minimize weight and maximize tube life.
- Affordable design for smaller budgets.



Cool!

Optional Pediatric Skins – Yellow Submarine, Friendly Jungle and Outer Space designs to help put your little patients at ease.

Meet the Revolution Nano.

Carestream's DRX-Revolution Mobile X-ray System changed the way the industry thinks about mobile imaging. Expanding the Carestream portfolio, we have now developed a fully integrated, digital, budget-friendly mobile system. The new DRX-Revolution Nano Mobile X-ray System is designed to deliver exceptional value and performance. The Revolution Nano utilizes Carbon Nano Tube Technology. The configuration of a Carbon Nano Tube (CNT) field-emitter tube differs significantly from traditional thermionic tubes. The CNT emitter in an X-ray tube consists of a large number of CNTs arranged vertically on a conductive substrate and replaces the traditional filament.

This technology allows the Revolution Nano to be smaller and more nimble than the DRX-Revolution, making it ideal to maneuver through small patient rooms, crowded hallways or tight spaces in the ICU or NICU.

The DRX-Revolution Nano is ideal for portable chest, ICU, orthopaedic and pediatric imaging. It delivers all of these key advantages:

- A unique Carbon Nano Tube features an ultra-lightweight design for easy maneuverability and arm positioning.
- Integrated DR technology delivers faster workflow and higher image quality than CR systems – at a lower capital cost.
- *Feather Touch Drive* and small footprint make maneuvering and navigating simple and easy.
- Sleek design offers line-of-sight visibility over the column and around the entire system.
- X-factor design allows detector sharing between rooms or mobile systems.
- Uses Eclipse image processing to power our software applications including Tube and Line Visualization, Pneumothorax Visualization, Bone Suppression, SmartGrid and Pediatric applications.

Everybody Benefits.

Everyone in your facility will appreciate the benefits of the DRX-Revolution Nano.

Radiologist

- High-quality diagnostic images: reliably and consistently
- Advanced software options for specialized visualizations

Radiology Administrator

- Lower capital investment
- Shareable detector to maximize productivity
- Future-proof design to prevent technology obsolescence
- Maximized return on investment

Radiographer

- Lightweight, easy to maneuver and position
- Clear line of sight
- Convenient bin design to put all supplies within reach
- Same intuitive user interface as all Carestream systems

Patients

- A quicker, more comfortable imaging process – no need to transport patients to X-ray rooms
- High-quality images to support accurate and timely diagnoses and treatment

A Fully Integrated DR Solution.

Driven by powerful CARESTREAM ImageView Software, the DRX-Revolution Nano can help speed your workflow, boost your productivity and improve your standard of care – all at an exceptionally affordable price.



Carbon Nano Tube with ultra-lightweight design and long five-year tube life

Balanced Articulating Arm – suitable for any kind of projection – unfolds and rotates for easy positioning

Lithium Iron-Phosphate Battery – longer lasting and quicker recharging than lead-acid batteries

Battery-Powered Exposure eliminates the need to plug in while capturing images

Ergonomically Designed Bin for storing supplies within easy reach

SmartGrid image-processing software, an image-enhancement algorithm to address the challenges caused by scatter and provide automatic processing to create reduced-scatter images



LED/Collimator Controls on tusks to maximize the radiographer's reach for greatest-possible SID

Feather Touch Drive for easy maneuverability

Independent Diaphragm Controls – no need to move the tube head

Detector-Bagging Feature makes bagging easy and helps reduce the possibility of detector drops

Kickstand to easily move over thresholds and elevators



Industry-Leading Performance: Meet Our Family of DRX Detectors.

The DRX-Revolution Nano is compatible with DRX-1, DRX Plus, DRX Core and DRX Plus 2530C Detectors.

NO MATTER WHICH ROAD YOU'RE ON,
WE'LL GET YOU WHERE
YOU NEED TO GO.

SMART IMAGING SOLUTIONS



The Fastest Route to a Successful Imaging Future.

As you travel to the future of imaging, what's your next step? Perhaps you need to accelerate your X-ray workflow. Or add a new imaging solution to your facility to provide better patient care.

Wherever you're headed, we have the smart solutions you need: Full-digital imaging rooms and leading-edge mobile units. Wireless, shareable detectors and CR-to-DR retrofit kits. Affordable CR systems. Advanced digital imagers and our Managed Print Solutions service. And as always, the gold-standard in X-ray film technology.

Start Mapping Your Route Today.

From small clinics, urgent care facilities and specialty practices to large hospitals and imaging centers, every facility will find a Carestream solution to meet the most specific needs. Visit carestream.com and start your journey to greater productivity and a higher standard of patient care.

Carestream also offers Non-Destructive Testing solutions to help ensure quality and safety in the aerospace, petrochemical and other industries, as well as Contract Manufacturing to bring innovative coated products to market.

A Community of Service and Support.

For dependable service, look to our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver.



Eclipse: Harness the Power of Decades of Innovation.

Carestream's Eclipse is the engine behind our innovative imaging software. It uses AI technology and proprietary algorithms to significantly increase the value of the entire imaging chain, from capture to diagnosis.

carestream.com/drx-nano



"Rx only"

© Carestream Health, Inc., 2021. CARESTREAM is a trademark of Carestream Health. CAT 2000 216 02/21

Carestream